

AACR Virtual Special Conference

TUMOR IMMUNOLOGY AND IMMUNOTHERAPY

October 5-6, 2021

AACR

American Association
for Cancer Research®

FINDING CURES TOGETHER®

**A promising cancer immunotherapy target:
Novel agonistic human antibody against the human
T-cell costimulatory receptor CD27**

Thierry Guillaudeux, PhD
Senior Vice President Immuno-Oncology
Kineta, Inc. Seattle, WA

tguillaudeux@kineta.us

A promising cancer immunotherapy target: Novel agonistic human antibodies against the human T-cell costimulatory receptor CD27



Yulia Ovechkina, Shaarwari Sridhar, Emily Frazier, Nathan Eyde, Neda Kabi, Jessica Cross, Carla Talbaux, Eric Tarcha, Shawn Iadonato and Thierry Guillaumeux
Kineta Inc. 219 Terry Av. North, Seattle, WA 98109.



Background

- CD27 is a member of the TNF-Receptor superfamily and plays a critical role in T-cell activation by providing a costimulatory signal after binding to its ligand CD70.
- CD27 signaling enhances T-cell proliferation, activation and differentiation of effector and memory T cells and therefore promotes cytotoxic T cell (CTL)-based anti-tumor immunity.
- CD27 is expressed on CD4+ and CD8+ T cells, B cells and NK cells.
- Agonistic stimulation of CD27 is a promising cancer immunotherapy approach to reinvigorate specific T cell driven anti-tumor responses.

Objectives

Evaluate 147 fully human monoclonal anti-CD27 antibodies (mAb)


- Measure antibodies' affinity to CD27.
- Test anti-CD27 antibodies binding to human, murine and primate CD27 to evaluate species cross-reactivity.
- Perform competition assay to test if soluble CD70-ECD disrupts anti-CD27 antibody binding.
- Assess antibody agonist activity using NFkB-Luc2/CD27 Jurkat reporter assay.
- Characterize the CD27-mediated costimulatory effects of anti-CD27 antibodies on human peripheral blood T cells.

Exceptional antibody diversity in both heavy and light chains


147 fully human ScFv antibodies directed against human CD27 obtained after immunization of Trianni mice and single B cell sequencing

- > 14 V_H diversity groups
- > 8 V_L diversity groups
- Highest diversity in CDR3-H

CDR3_Heavy Chains



CDR3_Light Chains



Kineta's anti-CD27 mAbs bind to human-CD27 with different affinity




Figure 1: (A) ELISA binding studies were conducted with immobilized human CD27 (rabbit-Fc). Binding of anti-CD27 mAb was detected using biotinylated anti-human IgG (light chain) followed by streptavidin HRP. **(B)** Octet kinetics studies were performed using 1x kinetics buffer (pH 7.2) and 10 µg/ml (67 nM) of anti-CD27 mAb loaded onto an Anti-human IgG Fc Capture (AHC) biosensor for 120 seconds. The loaded biosensor was incubated with 50nM of human CD27-ECD (His-Tag) for 240 seconds for the association step and then transferred into the kinetic buffer for 360 seconds to measure the dissociation. 1:1 curve fitting analysis was performed to determine k_a, k_d, and K_D.

Anti-CD27 mAbs cross-react with NHP-CD27 but not with mouse-CD27




Figure 2: Binding of anti-CD27 mAb to immobilized human CD27 (rabbit-Fc), non-human primate (NHP) CD27 (human-Fc), or mouse CD27 (mouse-Fc). Bound anti-CD27 mAbs were detected using biotinylated anti-human IgG (light chain) followed by streptavidin HRP.

Anti-CD27 mAbs prevent CD70 binding




Figure 3: Disruption of soluble human CD70-ECD (mouse Fc) binding to immobilized human CD27 (rabbit-Fc) by anti-CD27 mAb in a competition ELISA. CD70-ECD (mouse-Fc) at 0.11nM was spiked with anti-CD27 mAb (0.03-66.67nM) and incubated on the plate coated with CD27 for 1 hour. Bound CD70 was detected using peroxidase anti-mouse IgG (H+L). EC₅₀ and % Blockage of CD27/CD70 binding was measured.

Results

Selected anti-CD27 mAbs induce strong agonistic activity




Figure 4: Anti-CD27 agonist assay was carried out with NFkB-Luc2/CD27 Jurkat cells seeded at 2e5 cells/well in X-VIVO-15 medium in the 96-well plates pre-coated with 5 µg/mL anti-CD3 mAbs. Cells were then incubated with anti-CD27 antibodies or CD27 ligand (CD70-ECD) and incubated for 4 hrs at 37°C, 5% CO₂ humidified incubator. NFkB reporter induction was measured by luminescence in the presence of Bio-Glo reagent.

Selected anti-CD27 mAbs demonstrate robust T cell activation




Figure 5: T-cells from two healthy donors were stained with cell trace violet (CTV) and stimulated with 10 µg/mL of anti-CD27 antibodies and T-cell activator (2 µg/mL anti-CD3 and 5 µg/mL anti-CD28). The frequency of CTV labeled cells and cytokine release were examined after 4 days using flow cytometry and Luminex assay, respectively. **(A)** Representative histogram plots comparing T-cell proliferation for 2 selected anti-CD27 mAbs to T-cell stimulator without anti-CD27 antibody. **(B)** Proliferation data were normalized against the results obtained with T-cell activator without anti-CD27 antibody for each donor. **(C)** Cytokine data are represented as pg/mL.

Conclusions

Identified lead candidates for future pre-clinical development

- 147 fully human anti-CD27 monoclonal antibodies with unique sequences were generated.
- Anti-CD27 agonist assay showed strong agonist activity for 8 pre-selected anti-CD27 antibodies.
- Human T cell activation assay showed increased proliferation and cytokine secretion for only 4 mAbs.
- Further *in vitro* and *in vivo* developments are on-going to select our lead anti-CD27 agonist antibody.